



Preparation for 3D cell culture on Alvetex® Scaffold

1. CHO-K1 cells (ATCC, CCL-61) were routinely maintained in T-75 flasks.





Figure 1. Phase contrast micrographs of CHO-K1 cells grown in conventional 2D culture plates. Images show cells at low (left) and high (right) confluency. Scale bars: 100 μ m.

2. Complete medium consisted of: F-12K nutrient mixture (Kaighn's modification; Gibco 21127022) supplemented with 10 % v/v FBS and 100 μ g/ml Penicillin and 10 μ g/ml Streptomycin. *Note: Medium already contains L-glutamine.*
3. Cells were harvested by trypsinisation and centrifuged for 5 minutes (1000 rpm). The supernatant was discarded and the cell pellet was re-suspended in appropriate volume of media for cell counting by Trypan Blue.
4. Cells were re-suspended at a concentration of 1×10^6 cells / 150 μ l for seeding per well.
5. Alvetex® Scaffold 12-well discs were prepared for seeding by washing in 70 % ethanol and rinsing once with 4 ml of medium per well.
6. 150 μ l of the cell suspension was added to the centre of the Alvetex® Scaffold disc, which was equivalent to 1×10^6 cells per well.
7. The plate was incubated 30 minutes at 37 °C with 5 % CO₂ to allow the cells to settle into the scaffold.
8. 4 ml of medium was added to each well taking care not to dislodge cells from Alvetex® Scaffold.
9. Plates were re-incubated and maintained by complete medium exchange after every 1-2 days.



UK & Rest of World

184 Milton Park, Abingdon
OX14 4SE, Oxon, UK
Tel: +44 (0) 1235 828 200
Fax: +44 (0) 1235 820 482

Switzerland

Centro Nord-Sud 2E
CH-6934 Bioggio-Lugano
Tel: +41 (0) 91 604 55 22
Fax: +41 (0) 91 605 17 85

Deutschland

Bockenheimer Landstr. 17/19
60325 Frankfurt/Main
Tel: +49 (0) 69 779099
Fax: +49 (0) 69 13376880

United States

23591 El Toro Rd, Suite #167
Lake Forest, CA 92630
Tel: + 1 800 987 0985
Fax: + 1 949 265 7703




Figure 2. Brightfield micrographs showing the structure of CHO-K1 cells cultured for 3 days on 22 mm diameter Alvetex® Scaffold discs presented in 12-well plates. Cells were fixed, embedded in paraffin wax, sectioned (10 µm) and counterstained with Haematoxylin and Eosin.





Figure 3. Biochemical analysis of cell viability using a standard MTT assay. Data from 3 sample replicates of CHO-K1 cells are shown. Each well was sampled in triplicate. Mean values are shown. Cells were cultured for 3 days on 22 mm Alvetex® Scaffold discs presented in 12-well plates.