

Preparation for 3D cell culture on alvetex®

1. HepG2 cells were routinely maintained in T-75 flasks.




Figure 1. Phase contrast micrographs of HepG2 cells grown in conventional 2D culture plates. Images show cells at low (left) and high (right) confluency. Scale bars: 100 μ m.

2. Complete media consisted of: MEM media (Gibco 21090) supplemented with 10 % v/v FBS, 2 mM L-glutamine and 100 U/ml Penicillin/Streptomycin.
3. Cells were harvested by trypsinisation and centrifuged for 5 minutes (1000 rpm). The supernatant was discarded and the cell pellet was re-suspended in an appropriate volume of media for cell counting by Trypan Blue.
4. Cells were re-suspended at a concentration of 1.6×10^7 cells/ml for seeding.
5. Alvetex® 12-well plates were prepared for seeding with 70% ethanol (2 ml per well) and media washes (twice with 3 ml of media each) as described in the product information leaflet.
6. 125 μ l of the cell suspension was added to the centre of the alvetex® disc, which was equivalent to 2×10^6 cells per well.
7. The plate was incubated 3 hours at 37 °C with 5% CO₂ to allow the cells to settle into the scaffold.
8. 4 ml of media was added to each well taking care not to dislodge cells from alvetex®.
9. Plates were re-incubated and maintained by complete media exchange after every 1-2 days.






Figure 2. Brightfield micrographs showing the structure of HepG2 cells cultured for 7 days on 22 mm diameter alvetex® discs presented in the 12-well plate format. Cells were fixed, embedded in paraffin wax, sectioned (10 μm) and counterstained with haematoxylin and eosin.





Figure 3. Biochemical analysis of cell viability using a standard MTT assay. Data from 3 sample replicates of HepG2 cells are shown ($n=3$, mean \pm SD). Cells were cultured for 3 days on 22 mm alvetex® discs presented in the 12-well plate format.