




presented in different formats


Comparison of 3D cell growth patterns on alvetex[®] 12-well plate format (AMS.AVP002)

Alvetex® in the 12-well plate format (AMS.AVP002) was treated with EtOH and washed

complete medium prior to cell seeding. [Complete medium consisted of: DMEM, 10% FBS, 2 mM L-glutamine and 100 U/ml Penicillin & Streptomycin]. HaCaT cells (a human keratinocyte cell line) were plated at a density of 0.5×10^6 cells in 150 μ l per well while HepG2 cells (a human liver cell line) were seeded at 2x10⁶ cells in 150 µl per well. Plates were incubated for three hours before flooding with further media and maintained for 7 days. After preserving in Bouins fixative the discs were paraffin embedded, sectioned (10 µm) and counterstained with Haematoxylin and Eosin. HaCaT cultures demonstrated significant cell invasion into the matrix, while HepG2 cells remained resident in the top 25% of the matrix.



HaCaT cells



HepG2 cells

Figure 1. Comparison of 3D cell growth patterns on alvetex® 12-well plate format (AMS.AVP002). Micrographs taken at 20x magnification.




amsbio




Comparison of 3D cell growth patterns of HaCaT cells on alvetex® presented in various formats

HaCaT cells (a human keratinocyte cell line) were seeded (0.5×10^6 cells in 150 μ l per well) on EtOH-treated and complete medium washed alvetex[®] scaffolds in the following formats: 12-well plate (AMS.AVP002), 6-well inserts (AMS.AVP004-32) in 6-well plate and 12-well inserts (AMS.AVP005-34) in well insert holder in deep Petri dish (AMS.AVP015-2). Cultures were maintained for 7 days. [Complete medium consisted of: DMEM, 10% FBS, 2 mM L-glutamine and 100 U/ml Penicillin & Streptomycin]. After preserving in Bouins fixative the discs were paraffin embedded, sectioned (10 μ m) and counterstained with Haematoxylin and Eosin.



12 well plate (AMS.AVP002)




6-well insert (AMS.AVP004-32) in a 6-well plate







amsbio



12-well inserts (AMS.AVP005-34) in well insert holder in deep Petri dish (AMS.AVP015-2)

Figure 2. Comparison of 3D cell growth patterns of HaCaT cells grown on various alvetex® formats. Micrographs taken at 20x magnification.

Note significantly more proliferation and cell invasion into the alvetex[®] scaffold in cultures grown in well-inserts, (media interconnected feeding regime) due to nutrient availability from above and below the scaffold. Thus, as media volume and availability increased so did cell proliferation and scaffold penetration. In the case of well inserts contained in a well insert holder in a deep Petri dish, this resulted in the formation of a slab of tissue-like material.

