


CATALOG # 100332-1A
Package Size: 2 units / vial
Chondroitinase ABC Protease Free (*Proteus vulgaris*)
EC 4.2.2.4 / CAS Number: 9024-13-9
DESCRIPTION:

Chondroitin ABC lyase, Chondroitinase ABC, Chondroitin ABC eliminase.

Chondroitinase ABC Protease Free, catalyzes the removal of chondroitin sulfate and dermatan sulfate side chains of proteoglycans, yielding protein enriched core molecules with the enzymatically modified linkage oligosaccharides^{1,2}. The enzyme is highly specific for the galactosaminoglycan chains and without activity on core proteins, keratan sulfate chains, and heparin/heparan sulfate chains even in the absence of inhibitors for proteases, keratanases, and heparitinases.

Since the purity of the enzyme preparation is extremely high, and since no protein stabilizer is added, the enzyme preparation gives a main protein band with relative molecular weight of about 80,000 on SDS polyacrylamide gel electrophoresis. Because of these properties, Chondroitinase ABC Protease Free is suited for analysis and/or preparation of proteoglycan core proteins (for details of the procedures, see Ref. 1 and 2).

REACTION:


Where:	GalNAc	<i>S</i> - <i>N</i> -acetylgalactosamine 4 (or 6)-sulfate
	Hex A	D-glucuronic acid or L-iduronic acid
	Glc A	D-glucuronic acid
	Gal	D-galactose
	Xyl	D-xylose
	Ser	Serine
	ΔGlc A	Δ4,5-glucuronic acid
	ΔDi-S	unsaturated sulphated disaccharides

SPECIFICATIONS:

Activity:	≥2.0 units/vial	
Specific Activity:	≥110 units/mg protein	
Contaminants:	Chondro-4-sulfatase	<10 ⁻⁶ units/vial
	Chondro-6-sulfatase	<10 ⁻⁶ units/vial
	(By Morgan-Elson reaction)	
	Keratanase	Not detected
	Heparinase	Not detected
	Heparitinase	Not detected
	(0.4 unites of Chondroitinase ABC Protease Free, 37°C, 20 hours)	
	Protease	Not detected
	(By FITC-Casein ⁴)	

AMSBIO | www.amsbio.com | info@amsbio.com

 **UK & Rest of the World**
 184 Park Drive, Milton Park
 Abingdon OX14 4SE, UK
 T: +44 (0)1235 828 200
 F: +44 (0) 1235 820 482

 **North America**
 1035 Cambridge Street,
 Cambridge, MA 02141
 T: +1 (617) 945-5033 or
 T: +1 (800) 987-0985
 F: +1 (617) 945-8218

 **Germany**
 Bockenheimer Landstr. 17/19
 60325 Frankfurt/Main
 T: +49 (0) 69 779099
 F: +49 (0) 69 13376880

 **Switzerland**
 Centro Nord-Sud 2E
 CH-6934 Bioggio-Lugano
 T: +41(0) 91 604 55 22
 F: +41(0) 91 605 17 85

Appearance:	Lyophilised powder containing 20 mM Tris-HCl buffer, pH 7.2
Stabilizer:	BSA free
Preservative:	None
Reconstitution:	Dissolve the enzyme in 200 µl of 0.1% BSA
Optimum pH:	8.0 (for chondroitin sulfate) 6.2 (for hyaluronic acid)
Recommended Reaction Temperature:	37°C
Molecular Weight:	120,000-145,000 (gel filtration)
Activators:	Acetate (0.05 M)
Inhibitors:	Zn ²⁺ (10 ⁻³ M ZnCl ₂ inhibits by 100%) Heparin (an equimolar amount of heparin inhibits by ca. 70%)

Unit Definition:

One unit is defined as the quantity of the enzyme that catalyzes the formation of 1 µmole of unsaturated disaccharide from chondroitin sulfate C (shark cartilage, see product code **400670** on the website) per minute at 37°C, pH 8.0.

ASSAY FOR ENZYME ACTIVITY:

Method:	The assay is based on that of Yamagata et al. ³	
Reaction mixture		
Substrate and Buffer solution:	0.1 µmole chondroitin sulfate C (shark cartilage, Cat. No. 400670 in:	
	0.4 M Tris-HCl buffer, pH 8.0	10 µl
	0.4 M sodium acetate	10 µl
	0.1% bovine serum albumin (BSA)	10 µl
	Distilled water	70 µl
Enzyme solution:	Diluted enzyme (1-5 mU) with 0.1% BSA	20 µl
Total volume:		120 µl

Procedure

Reaction: The reaction mixture is incubated at 37°C for 10 minutes and stopped by boiling for 1 minute.

Morgan-Elson Reaction: To enzyme reaction mixture (120 µl) add 0.1 ml of 5% K₂B₄O₇ (pH 9.0), and heat in a boiling water bath for 7 minutes. After cooling, add 1 ml of glacial acetic acid, mix, add 0.4 ml of the Morgan-Elson reagent, and incubate at 37°C for 20 minutes. Measure A₅₈₅.

Calculation:

$$\text{Enzyme Unit (units/ml)} = \frac{A_{585}}{2.38} \times \frac{0.1}{G} \times \frac{1}{10} \times \frac{1}{E}$$

where: G: Adsorption of 0.1 µmole GalNAc (A₅₈₅)
E: Volume of enzyme solution (ml)

Note:

Since the enzyme action is inhibited more or less by keratan sulfate, dermatan sulfate, or heparin/heparin sulfate, amounts of the enzyme to be used should be determined by pilot tests.

STORAGE:

Store at -70°C until opened. Following reconstitution, aliquot and freeze at -20°C.

REFERENCES:

- 1) Oike, Y., Kimata, K., Shinomura, T., Nakazawa, K. and Suzuki, S. (1980) *Biochem. J.*, **191**, 193
- 2) Oike, Y., Kimata, K., Shinomura, T., Suzuki, S., Takahashi, N. and Tanabe, K. (1982) *J. Biol. Chem.*, **257**, 9751
- 3) Yamagata, T., Saito, H., Habuchi, O. And Suzuki, S. (1968) *J. Biol. Chem.*, **243**, 1523
- 4) Twining, S.S. (1984) *Anal. Biochem.*, **143**, 30

NOTE:

For *in vitro* research use only – not for diagnostic or therapeutic use. This product is not a medical device.

AMSBIO | www.amsbio.com | info@amsbio.com

 **UK & Rest of the World**
184 Park Drive, Milton Park
Abingdon OX14 4SE, UK
T: +44 (0)1235 828 200
F: +44 (0) 1235 820 482

 **North America**
1035 Cambridge Street,
Cambridge, MA 02141
T: +1 (617) 945-5033 or
T: +1 (800) 987-0985
F: +1 (617) 945-8218

 **Germany**
Bockenheimer Landstr. 17/19
60325 Frankfurt/Main
T: +49 (0) 69 779099
F: +49 (0) 69 13376880

 **Switzerland**
Centro Nord-Sud 2E
CH-6934 Bioggio-Lugano
T: +41(0) 91 604 55 22
F: +41(0) 91 605 17 85