

# Instruction Manual: StemFit® Basic03 - Clinical Grade

## For maintenance and expansion of human ES/iPS cells



### 1. Materials Required

StemFit® Basic03 - Clinical Grade (amsbio #SFB-503)

Cell dissociation reagents (e.g. Accutase, Detachin #T100100),

Extracellular Matrix (ECM) - hESC qualified (Recommended iMatrix-511 amsbio #AMS.892)

Human bFGF (amsbio #AMS-480)

Y-27632

PBS (-)

### 2. Media Preparation

StemFit® Basic03 (Basic03) is provided frozen as a 2-component set containing “Liquid A” and “Liquid B”, and can be stored at below -20°C until use. Use sterile techniques to prepare Basic02 medium.

- 1) Before use, thaw frozen “Liquid A” and “Liquid B” with occasional mixing at room temperature (15-25°C).  
**CAUTION: Do not thaw “Liquid B” at 37°C,** as it accelerates the degradation of the medium ingredients.
- 2) Aseptically mix medium components by adding the full volume of “Liquid B” to “Liquid A”. Mix thoroughly.
- 3) Upon thawing, StemFit® Basic03 medium can be aseptically aliquoted and stored at below -20°C. Before use, thaw an aliquot in the refrigerator overnight.

- 4) **Add bFGF at a concentration of 10 ng/ml.**

**Note:** We recommend adjusting the concentration of bFGF (e.g. 40 - 80 ng/ml) according to suit your cell line if your cells differentiate.

- 5) Store the thawed medium in the refrigerator.

**CAUTION: Thawed StemFit® Basic03 medium may be stored at 2 - 8 °C for up to two weeks.**

**CAUTION: We recommend storing the medium in the dark.**

- 6) Before use, warm aliquots to room temperature and use immediately.

**CAUTION: Do not heat the thawed medium to 37°C.**

### 3. Passage Protocol (6-well plate; Also see our technical tips: Key points for successful single-cell passage)

- 1) Culture vessel coating: Add LDEV-free hESC-qualified ECM to cold DMEM/F-12 at a 1:100 ratio and mix well immediately. Add 1 ml of the ECM mixture to one well of a six-well plate. Incubate at 37°C for at least 1 h.

**Note:** You can use other matrices such as iMatrix-511 laminin-511, Matrigel, vitronectin, laminin-521 or laminin-511E8.

- 2) Cell passage: Aspirate the medium and wash once with 2 ml of PBS/well.
- 3) Aspirate the PBS and add 500 µl of Accutase. Incubate at 37°C for 10 min.

**Note:** TrypLE can also be used for cell dissociation.

**Note:** Incubation time may vary depending on the matrix.

- 4) Pipette the cells to fully dissociate and transfer cells to a 15-ml tube filled with 500 µl of Basic02 supplemented with bFGF (Basic02+F) containing Y-27632 (final concentration: 10 µM).
- 5) Count the cells with a cell counter or hemocytometer (optimized for the cell types).
- 6) Centrifuge the tubes at 300 g at room temperature for 4 min.


- 7) Aspirate the medium and resuspend cells at a density of 1,000 cells/ $\mu$ l.
- 8) Aspirate the Geltrex solution and add 1.5 ml of Basic02+F containing Y-27632/well (final concentration: 10  $\mu$ M).
- 9) Add 10-20  $\mu$ l of resuspended cells directly to the new well (10,000-20,000 cells/well).
- 10) Culture the cells at 37°C in a 5%CO<sub>2</sub> incubator >24 hours.
- 11) Aspirate the medium and add 1.5 ml of Basic02+F
- 12) Perform medium changes with 1.5 ml of Basic02+F.
- 13) Passage the cells every 7 d.

Note: You can culture hPSCs without weekend medium changing. See the following passage schedule examples.

**CAUTION: If the color of the medium turns orange or yellow, it should be changed every day.**

**CAUTION: Do not allow cells to become confluent.**


#### StemFit® Passage Schedule Examples (Weekend-free feeding)



#### 4. Transfer from other culture systems

- To transfer cells from other culture systems to the StemFit® system, we recommend passaging with the original culture system then switching the culture medium to Basic02 supplemented with bFGF (Basic02+F) 2 - 3 days prior to the next passage.
- Seeding the cells at a higher density (>1.0 x 10<sup>5</sup> cells per well (6-well plate)) may be helpful for the first few passages.

#### Transition Schedule Example



#### 5. Reference

Morizane, R. & Bonventre, J. V. Generation of nephron progenitor cells and kidney organoids from human pluripotent stem cells. *Nat Protoc.* 2017 Jan;12(1):195-207.

## 6. FAQs & Troubleshooting

- 1) What are the benefits of single cell culture? / Why is single cell culture recommended?
  - High fold expansion rate (~100x expansion / weekly passage)
  - Reproducible and manageable culture by controlling the numbers of seeded cells
  - Cost-effective culture with lower medium volume and less frequent medium changes
  - Produce an iPSC colony derived from single cells. (essential for genome editing)
  
- 2) Can I use StemFit® for clump culture?
  - Yes, but we recommend making a small clump and seeding at a low cell density.
  
- 3) Cells do not grow well.
  - Adjust the bFGF concentration (e.g. 40 - 80 ng/ml) according to your cell line
  - Try a higher seeding density (e.g. > 1.0 x 10<sup>5</sup> cells per well (6-well plate))
  - Distribute the cells evenly upon passage
  - Culture in Y-27632-containing medium for more than 24 hours
  - Make sure that the medium was thawed within 2 weeks and has not been heated to 37°C

AMSBIO | [www.amsbio.com](http://www.amsbio.com) | [info@amsbio.com](mailto:info@amsbio.com)

 **UK & Rest of the World**  
184 Park Drive, Milton Park  
Abingdon OX14 4SE, UK  
T: +44 (0)1235 828 200  
F: +44 (0) 1235 820 482

 **North America**  
1035 Cambridge Street,  
Cambridge, MA 02141  
T: +1 (617) 945-5033 or  
T: +1 (800) 987-0985  
F: +1 (617) 945-8218

 **Germany**  
Bockenheimer Landstr. 17/19  
60325 Frankfurt/Main  
T: +49 (0) 69 779099  
F: +49 (0) 69 13376880

 **Switzerland**  
Centro Nord-Sud 2E  
CH-6934 Bioggio-Lugano  
T: +41(0) 91 604 55 22  
F: +41(0) 91 605 17 85