

I. Introduction

StemFit® Basic04 culture medium.

NOTE: EXPERIMENT SCALE

In this protocol, instructions are described for a 35-mm dish or one well of a 6-well plate. For other well sizes, multiply the protocol's recipe by the ratio between well sizes.

II. Materials Required

- StemFit® Basic04 medium (amsbio, Cat#SFB-504)
- iMatrix-511 (Cat# AMS.892)
- bFGF (i.e. FGF2) growth factor (10 µg/ml, sterile-filtered) (amsbio, Cat#AMS-480)
- ROCK inhibitor
- TrypLE Select (Life Tech, Cat#12563-029)
- Phosphate-buffered saline (PBS without Ca⁺⁺ Mg⁺⁺)
- 0.5 mM EDTA in PBS
- STEM-CELLBANKER (amsbio, Cat#11890 and Cat#11897)

III. Pre-Protocol Preparation

Complete StemFit® Medium

1. Thaw StemFit® at 4°C overnight. Before use, thaw at room temperature (15-25°C). CAUTION: Do not thaw at 37°C, as this accelerates the degradation of the medium components.
2. Upon thawing, add bFGF at a concentration of 10 ng/ml.

Note: For new users and cell lines we recommend starting with a concentration of bFGF (e.g. 40 - 80 ng/ml) according to suit your cell line if your cells differentiate.

CAUTION: Thawed StemFit® Basic04 medium may be stored at 2 - 8°C for up to two weeks.

CAUTION: We recommend storing the medium in the dark

3. StemFit® Basic04 medium can be aseptically into 45 ml aliquots in 50 ml Conical Centrifuge Tubes (e.g., Falcon™, Corning®), store aliquots at -20°C. Before use, thaw an aliquot in the refrigerator overnight. Warm aliquots to room temperature and use immediately.

CAUTION: Do not heat the thawed medium to 37°C.

Coated Plates with ECM

NOTE: CHOICE OF MATRICE

For maximum cloning efficiency, we recommend using recombinant laminin-511 E8 available from amsbio. StemFit® Basic04 has been proven to perform with alternative matrices such as Matrigel, vitronectin, Geltrex and laminin-521

1. Take 12 ml ice-cold PBS into a tube and add 28.5 µl iMatrix-511 (0.5 mg/ml). Mix them well. This mixture is stable for up to 2 weeks at 4°C.

2. Add 2 ml diluted iMatrix to one well of a 6-well plate or a 35-mm culture dish.
3. Incubate the plate or dish at 37°C, 5% CO₂ for 2 hours (or 4°C overnight.)
4. Aspirate the supernatant from each well or dish and add 2 ml PBS to it.
5. Incubate the plate or dish at 37°C, 5% CO₂ until user's hPSCs are ready for plating or up to 1 hour.

10 mM ROCK Inhibitor

1. Prepare ROCK inhibitor stocks at a concentration of 10 mM in molecular-biology-grade DMSO. The stock is stable for up to 2 weeks at 4°C.
2. We recommend aliquoting 10 mM ROCK Inhibitor into 100 µl aliquots and storing the aliquots at -20°C.

0.5x TrypLE Select

1. Dilute TrypLE-Select (1x) with 0.5 mM EDTA in PBS at a ratio of 1:1. The mixture is referred to as 0.5x TrypLE Select and is stable for up to 2 months at 4°C.

Note: A trypsin inhibitor is not required when using TrypLE Select as dilution will inactivate TrypLE Select mostly. However, we still recommend minimizing carryovers of TrypLE Select in culture.

IV. Protocol

NOTE: CELL CONFLUENCY

hPSC culture is best for passaging when the cell confluency reaches 80-90%. However, when the growing colonies start to touch each other, consider harvesting them even if the cell confluency is below 80%. Otherwise, cells begin to differentiate.

Day 1-2 - Passaging

1. Warm Complete StemFit® Medium and 0.5x TrypLE Select for 20-30 minutes to room temperature.
2. Take 4 ml Complete StemFit® Medium into a tube and add 4 µl 10 mM ROCK Inhibitor to it. Mix them well.
3. Aspirate old medium from the culture and add 1.5 ml Complete StemFit® Medium with 10 µM ROCK inhibitor to it.
4. Incubate the culture at 37°C, 5% CO₂ for 1 hour. This treatment enhances the survival of hPSCs.
5. Aspirate the old medium from the culture and add 2 ml PBS to it.
6. Aspirate PBS from the culture and add 300 µl 0.5x TrypLE Select to it.
7. Incubate the culture at 37°C, 5% CO₂ for up to 10 minutes.
8. Carefully aspirate 0.5x TrypLE Select from the culture using a P1000 pipettor and add 1 ml PBS to it.
9. Gently rock the plate or dish 3 times.
10. Carefully pipet out PBS using a P1000 pipettor and add 1 ml Complete StemFit® Medium with 10µM ROCK inhibitor.
11. Disperse the medium over the well bottom surface by pipetting 8-15 times to detach cells.
12. Collect the cell suspension in a tube.
13. Count cells to determine the volume of cell suspension needed to plate 1×10^4 cells.
14. Take out the determined volume of the cell suspension from the previous step in a tube and bring up the volume to 1.5 ml with StemFit® with 10 µM ROCK inhibitor.
15. Aspirate PBS from a new iMatrix-511 coated well and add 1.5 ml cell suspension to it.
16. Rest of the cells can be used for assays or frozen down using STEM-CELLBANKER (2×10^5 cells per vial).

17. Incubate the culture at 37°C, 5% CO₂ for 1 hour.
18. Observe the culture under a microscope and make sure that cells are evenly distributed in the well or dish.
19. Incubate the culture at 37°C, 5% CO₂ for 2 days.

NOTE: CONFLUENCY AND FEEDING FREQUENCY

Data indicates that a culture prepared using 1.0×10^4 cells reaches 80-90% confluency and produces about $2-3 \times 10^6$ cells in 7-8 days.

When 1.0×10^4 cells are plated on each well or dish, feed cultures every 2 days for the first 4 days.

When more than 1.0×10^4 cells are plated on each well or dish, feed cultures with 2-4 ml medium every day.

The culture reaches 80-90% confluency in 3-4 days.

Day 3-4 - Feeding

1. Warm Complete StemFit® Medium for 20-30 minutes to room temperature.
2. Aspirate the old medium from each well or dish and add 1.5 ml Complete StemFit® Medium to it.
3. Incubate the culture at 37°C, 5% CO₂ for 2 days.

Day 5-7 - Feeding

1. Warm Complete StemFit® Medium for 20-30 minutes to room temperature.
2. Aspirate the old medium from each well or dish and add 1.5 ml Complete StemFit® Medium to it.
3. Incubate the culture at 37°C, 5% CO₂ overnight.
4. Repeat steps 1-3 until Day 8.

CAUTION: If the color of the medium turns orange or yellow, it should be changed every day.

CAUTION: Do not allow cells to become confluent.

V. FAQs & Troubleshooting

- 1) What are the benefits of single cell culture? / Why is single cell culture recommended?
High fold expansion rate (~100x expansion / weekly passage)
Reproducible and manageable culture by controlling the numbers of seeded cells
Cost-effective culture with lower medium volume and less frequent medium changes
Produce an iPSC colony derived from single cells. (essential for genome editing)
- 2) Can I use StemFit® for clump culture?
Yes, but we recommend making a small clump and seeding at a low cell density.
- 3) Cells do not grow well.
Adjust the bFGF concentration (e.g. 40 - 80 ng/ml) according to your cell line
Try a higher seeding density (e.g. $> 1.0 \times 10^5$ cells per well (6-well plate))

Distribute the cells evenly upon passage

Culture in Y-27632-containing medium for more than 24 hours

Make sure that the medium was thawed within 2 weeks and has not been heated to 37°C

VI. References

1. Nakagawa M, Taniguchi Y, Senda S, Takizawa N, & Ichisaka T (2014). A novel efficient feeder-free system for the derivation of human induced pluripotent stem cells. *Scientific Reports*, 8: 1-7.
2. Akiyama T, Wakabayashi S, Soma A, Sato S, Nakatake Y, Oda M, Murakami M, Sakota M, Chikazawa-Nohtomi N, Ko SB, Ko MS (2016). Transient ectopic expression of the histone demethylase JMJD3 accelerates the differentiation of human pluripotent stem cells. *Development*, 143: 3674-3685.
3. Goparaju SK, Kohda K, Ibata K, Soma A, Nakatake Y, Akiyama T, Wakabayashi S, Matsushita M, Sakota M, Kimura H, Yuzaki M, Ko SB & Ko MS (2017). Rapid differentiation of human pluripotent stem cells into functional neurons by mRNAs encoding transcription factors. *Scientific Reports* 7, 42367.
4. Morizane R, Bonventre JV (2017). Generation of nephron progenitor cells and kidney organoids from human pluripotent stem cells. *Nat Protoc.*, 12: 195-207.

For Research Use Only. Not Intended For Human or Animal Diagnostic or Therapeutic Uses

Version 1 (May, 2019)

AMSBIO | www.amsbio.com | info@amsbio.com



UK & Rest of the World
184 Park Drive, Milton Park
Abingdon OX14 4SE, UK
T: +44 (0)1235 828 200
F: +44 (0) 1235 820 482



North America
1035 Cambridge Street,
Cambridge, MA 02141
T: +1 (617) 945-5033 or
T: +1 (800) 987-0985
F: +1 (617) 945-8218



Germany
Bockenheimer Landstr. 17/19
60325 Frankfurt/Main
T: +49 (0) 69 779099
F: +49 (0) 69 13376880



Switzerland
Centro Nord-Sud 2E
CH-6934 Bioggio-Lugano
T: +41(0) 91 604 55 22
F: +41(0) 91 605 17 85