



Embryonic Body Formation on LIPIDURE®-COAT




Materials

- Lipidure®-Coat Plate A-U96 (Cat No.AMS.51011610)
- mES cells of the cell line 129SV
- Medium for EB formation

IMDM or DMEM (Invitrogen), 20% FBS, 1 mM Sodium Pyruvate (Invitrogen), 0.1 mM non essential amino acids, 0.1 mM 2-mercaptoethanol, 50 U/mL penicillin, 50 mg/mL streptomycin



Method




- Step 1: ①Prepare of 80% confluent mES cells in 6 well plate.
 ②Wash by PBS at 3 times.
 ③Add 0.5~1 mL of 0.1% trypsin solution and incubate for 5 min at 37°C.
 ④Add 4 mL EB formation medium.
 ⑤Centrifuge (1,000 rpm, 5 min, 4°C).
 ⑥Discard supernatant and resuspend by EB formation medium.
 ⑦Count cell number and dilute into 5,000 cells/mL.

Step 2: ⑧Seed 200 µL/well (1,000 cells/well) in Lipidure®-Coat Plate A-U96.

Step 3: ⑨Incubate in the humidified CO₂-incubator (37°C, 5% CO₂) for 3~5 days.

Step 4: ⑩Harvest EB by pipetting.



UK & Rest of World

184 Milton Park, Abingdon
OX14 4SE, Oxon, UK
Tel: +44 (0) 1235 828 200
Fax: +44 (0) 1235 820 482

Switzerland

Centro Nord-Sud 2E
CH-6934 Bioggio-Lugano
Tel: +41 (0) 91 604 55 22
Fax: +41 (0) 91 605 17 85

Deutschland

Bockenheimer Landstr. 17/19
60325 Frankfurt/Main
Tel: +49 (0) 69 7790999
Fax: +49 (0) 69 13376880

North America

23591 El Toro Rd, Suite #180
Lake Forest, CA 92630
Tel: +1 800 987 0985
Fax: +1 949 265 7703

amsbio

info@amsbio.com

www.amsbio.com
AMS Biotechnology