

Spheroid Formation on LIPIDURE®-COAT

Materials

1. LIPIDURE®-COAT PLATE A-U96 (Cat No. AMS.51011610)
2. Cell; Rat Hippocampus Neuronal Cell (DS pharma Biomedical Cat No. MB-X0402D)
3. Human recombinant EGF
4. Human recombinant bFGF

Cell culture medium

1. Trial Set for Neuronal Cell Culture (DS pharma Biomedical Cat No. MB-X0601D)

Cell preparation

1. Neuronal cells were cultured according to Trial Set for Neuronal Cell Culture.

Spheroid formation


1. The cells were enzyme-treated and plated at a seeding density of 200, 1000, 5000 cells/well (0.1 mL/well) on LIPIDURE®-COAT Plate.
 2. The plate was incubated at 37°C in 5 %CO₂ incubator for 1~7 days.
 3. From day 1, the single spheroid were formed in each well.
- * For spheroid formation, the media consisting with 50ng/mL human recombinant EGF and 25ng/mL human recombinant bFGF was used.
- ** The medium was gently changed every other day by pipetting. (We recommend the change of 70 μ L/well.).

Application


1. You can use the spheroids for cell-based various assays.

Neurospheres formed on LIPIDURE®-COAT on day 5


Seeding density: 200 cells/well



Seeding density: 1000 cells/well



Seeding density: 5000 cells/well



Scale bar= 200 μm

All photos were supplied by Dr. Ijima at Kyusyu University



UK & Rest of World

184 Milton Park, Abingdon
OX14 4SE, Oxon, UK
Tel: +44 (0) 1235 828 200
Fax: +44 (0) 1235 820 482

Switzerland

Centro Nord-Sud 2E
CH-6934 Bioggio-Lugano
Tel: +41 (0) 91 604 55 22
Fax: +41 (0) 91 605 17 85

Deutschland

Bockenheimer Landstr. 17/19
60325 Frankfurt/Main
Tel: +49 (0) 69 779099
Fax: +49 (0) 69 13376880

North America

23591 El Toro Rd, Suite #180
Lake Forest, CA 92630
Tel: +1 800 987 0985
Fax: +1 949 265 7703

amsbio
info@amsbio.com
www.amsbio.com
AMS Biotechnology