

TREVIGEN[®] Product Data

For Research Use Only. Not For Use In Diagnostic Procedures.

amsbio

T4-Endo V FLARE[™] Module

Catalog #: 4055-100-FM

Contents:

Component	Catalog #	Amount	Storage
T4-Endo V	4055-100-01	100,000 units	4°C
25X FLARE [™] Buffer 1	3950-040-01	40 ml	Room Temp.
REC [™] Dilution Buffer	3950-010-03	10 ml	-20°C or 4°C
100X BSA Additive	3950-100-04	100 µl	-20°C or 4°C

Description: Reagents for analysis of UV DNA damage in single cells using the Trevigen CometAssay[™], FLARE Slides, and T4-Endo V.

Assay & Analysis: Refer to T4-Endo V FLARE Assay Kit's (Cat# 4055-100-FK) Instructions for Use for further details.

Note: CometSlides and FLARE Slides are interchangeable in the FLARE Assay. The difference is in the sample number; CometSlide holds 2 samples whereas the FLARE Slide holds 3.


Lot Specific Data: Contact Trevigen for lot specific activity of T4-Endo V.

© 2002 Trevigen, Inc. All Rights Reserved. Trevigen is a registered trademark and CometAssay, CometSlide, FLARE, and REC are trademarks of Trevigen, Inc. v21217

The product accompanying this document is intended for research use only and is not intended for diagnostic purposes or for use in humans.



Please
Recycle



UK & Rest of World
184 Milton Park, Abingdon
OX14 4SE, Oxon, UK
Tel: +44 (0) 1235 828 200
Fax: +44 (0) 1235 820 482

Switzerland
Centro Nord-Sud 2E
CH-6934 Bioggio-Lugano
Tel: +41 (0) 91 604 55 22
Fax: +41 (0) 91 605 17 85

Deutschland
Bockenheimer Landstr. 17/19
60325 Frankfurt/Main
Tel: +49 (0) 69 779099
Fax: +49 (0) 69 13376880

North America
23591 El Toro Rd, Suite #180
Lake Forest, CA 92630
Tel: + 1 800 987 0985
Fax: + 1 949 265 7703

amsbio
info@amsbio.com
www.amsbio.com
AMS Biotechnology