

# TREVIGEN® Product Data

For Research Use Only. Not For Use In Diagnostic Procedures



## UVDE FLARE™ Module

Catalog #: 4100-100-FM

### Contents:

Component	Catalog #	Amount	Storage
<i>S. pombe</i> UVDE GST-Δ228	4100-100-01	100 units	-80°C*
25X FLARE™ Buffer 2	3951-040-01	40 ml	Room temp.
REC™ Dilution Buffer	3950-010-03	10 ml	-20°C
100X BSA Additive	3950-100-04	100 µl	-20°C
100X Cation Solution	3950-100-05	100 µl	-20°C

\*For optimal stability, freeze in working aliquots at -80°C to avoid repeated freeze-thawings.

**Description:** Reagents for analysis of UV-induced DNA damage in single cells using the Trevigen CometAssay®, FLARE™ Slides, and *S. pombe* UVDE.

**Assay & Analysis:** Refer to the UVDE FLARE™ Assay Kit's (cat#: 4100-100-FK) Instructions for Use for further details.

Note: CometSlides and FLARE™ Slides are interchangeable in the FLARE™ Assay. The difference is in the sample number; the CometSlide™ holds 2 samples whereas the FLARE™ Slide holds 3.

**Lot Specific Data:** Contact Trevigen for lot specific activity of UVDE.

### Related Products:

#### FLARE™ Assay Kits:

Catalog #	Description	Damage Recognized	Size
4040-100-FK	<b>Fpg Kit</b>	8-oxoguanine, DNA containing formamidopyrimidine moieties	<b>75 samples</b>
4040-100-FM			<b>100 samples</b>
4055-100-FK	<b>T4 Endo V Kit</b>	Cis-syn isomers of cyclobutane pyrimidine dimers	<b>75 samples</b>
4055-100-FM			<b>100 samples</b>
4065-100-FK	<b>cvPDG Kit</b>	Cis-syn and trans-syn isomers of cyclobutane pyrimidine dimers	<b>75 samples</b>
4065-100-FM			<b>100 samples</b>
4130-100-FK	<b>hOGG1 Kit</b>	8-oxoguanine, DNA containing formamidopyrimidine moieties	<b>75 samples</b>
4130-100-FM			<b>100 samples</b>
4100-100-FK	<b>UVDE Kit</b>	Cyclobutane pyrimidine dimers, (6-4) photoproducts	<b>75 samples</b>
			<b>100 samples</b>
4045-010-FK	<b>Endonuclease III Kit</b>	Thymine Glycol, 5,6-dihydrothymine, urea, 5-hydroxy-6-hydrothymine, 5,6-dihydro-uracil, alloxan, 5-hydroxy-6-hydro-uracil, uracil glycol, 5-hydroxy-5-methylhydantoin, 5-hydroxycytosine, 5-hydroxy uracil, methyltartonyl urea, thymine ring saturated or fragmentation product	<b>75 samples</b>
4045-010-FM			<b>100 samples</b>

#### CometAssay® Kits:

Catalog #	Description	Size
4250-050-ESK	CometAssay® Starter Kit	each
4250-050-ES	CometAssay® Electrophoresis System	1 system
4250-050-K	CometAssay®	50 samples
4252-040-K	CometAssay® Higher Throughput Kit	40 samples
4254-200-K	CometAssay® Silver Staining Kit	200 samples
4253-096-K	CometAssay® Kit 96 Wells	96 samples

www.amsbio.com - info@amsbio.com

 <small>UK and rest of the world</small> Tel: +44 (0) 1235 828 200 Fax: +44 (0) 1235 820 482	 <small>DE</small> Tel: +49 (0) 69 779099 Fax: +49 (0) 69 13376880
 <small>CH</small> Tel: +41 (0) 91 604 55 22 Fax: +41 (0) 91 605 17 85	 <small>us</small> Tol Free: +1 800 987 0985 Tel: +1 949 265 7717 Fax: +1 949 265 7703

**amsbio**