
READY TO USE LENTIVIRUS PARTICLES

Price List

The logo for amsbio, featuring the word "amsbio" in white lowercase letters on a dark blue rectangular background. A small red dot is positioned above the letter "i".

amsbio

Fluorescent Markers

Code	Description	Pack Size	GBP	CHF	EUR
LVP001	gfp lentiviral / lentivirus particles	1x10E8 IFU/ml x 200ul	220	380	260
LVP001-PBS	GFP particles, in vivo ready	1x10E9 IFU/ml x 200ul	725	1260	860
LVP001-SF	Serum-free GFP particles	1x10E8 IFU/ml x 200ul	260	450	305
LVP002	gfp(his) lentiviral / lentivirus particles	1x10E8 IFU/ml x 200ul	220	380	260
LVP002-SF	Serum-free GFP(His) particles	1x10E8 IFU/ml x 200ul	260	450	305
LVP024	GFP inducible control particles	1x10E8 IFU/ml x 200ul	220	380	260
LVP011	cfp lentiviral / lentivirus particles	1x10E8 IFU/ml x 200ul	220	380	260
LVP011-SF	Serum-free CFP particles	1x10E8 IFU/ml x 200ul	260	450	305
LVP012	yfp lentiviral / lentivirus particles	1x10E8 IFU/ml x 200ul	220	380	260
LVP012-SF	Serum-free YFP particles	1x10E8 IFU/ml x 200ul	260	450	305
LVP023	rfp lentiviral / lentivirus particles	1x10E8 IFU/ml x 200ul	220	380	260
LVP023-PBS	rfp lentiviral / lentivirus particles, in vivo ready	1x10E9 IFU/ml x 200ul	725	1260	860
LVP023-SF	Serum-free RFP particles	1x10E8 IFU/ml x 200ul	260	450	305

Marker proteins & Enzymes

Code	Description	Pack Size	GBP	CHF	EUR
LVP009	luciferase lentiviral / lentivirus particles	1x10E7 IFU/ml x 200ul	220	380	260
LVP009-PBS	luciferase lentiviral / lentivirus particles, in vivo ready	5x108 IFU/ml x 200ul	725	1260	860
LVP009-SF	serum-free, luciferase lentiviral / lentivirus particles	5x10E7 IFU/ml x 200ul	260	450	305
LVP010	lacz lentiviral / lentivirus particles	1x10E7 IFU/ml x 200ul	220	380	260
LVP016	t-antigen lentiviral / lentivirus particles	1x10E7 IFU/ml x 200ul	220	380	260
LVP017-Bsd	tetr (bsd) lentiviral / lentivirus particles	1x10E7 IFU/ml x 200ul	220	380	260
LVP018	tk lentiviral / lentivirus particles	1x10E7 IFU/ml x 200ul	220	380	260
LVP019	tat lentiviral / lentivirus particles	5x10E7 IFU/ml x 200ul	220	380	260
LVP020	luciferase (his) lentiviral / lentivirus particles	1x10E7 IFU/ml x 200ul	220	380	260
LVP021	gfp_lacz (his) lentiviral / lentivirus particles	1x10E7 IFU/ml x 200ul	260	450	305
LVP021-SF	gfp_lacz (his) lentiviral / lentivirus particles	5x10E7 IFU/ml x 200ul	370	645	440
LVP013	cre lentiviral / lentivirus particles	1x10E7 IFU/ml x 200ul	220	380	260
LVP013-PBS	cre lentiviral / lentivirus particles, in vivo ready	5x108 IFU/ml x 200ul	725	1260	860
LVP013-SF	serum-free cre lentiviral / lentivirus particles	5x10E7 IFU/ml x 200ul	260	450	305
LVP014	flp recombinase lentiviral / lentivirus particles	1x10E7 IFU/ml x 200ul	220	380	260

Miscellaneous Lentiviral products

Code	Description	Pack Size	GBP	CHF	EUR
TLV-C	293t lentiviral / lentivirus packing cells	2 x 10E6 cell/ml x 1ml	220	380	260
LTSH-GB	pEco-Lenti-shRNA-(GFP-Bsd) cloning kit	1 Kit	480	835	570
LTSH-GP	pEco-Lenti-shRNA-(GFP-Puro) cloning kit	1 Kit	480	835	570
CC02	Eco(tm) _DH5a_Lenti chemical competent cells	20 x 50ul	85	145	100
HT-pack	lentiviral / lentivirus packaging plasmids	600ng/ul x 100ul	365	640	435

Human proteins for inducing pluripotent stem cells (iPSC)

Code	Description	Pack Size	GBP	CHF	EUR
LVP-stems	Full set (6 human stem factors) particles	200ul x 6	835	1455	990
LVP-stems-SF	Full set (6 human stem factors) serum-free particles	200ul x 6	1205	2100	1435
LVP003	hOCT3/4 inducible lentiviral / lentivirus particles	5x10E7 IFU/ml x 200ul	220	380	260
LVP003-SF	serum-free hOCT3/4 inducible particles	1x10E8 IFU/ml x 200ul	260	450	305
LVP004	hSOX2 inducible lentiviral / lentivirus particles	5x10E7 IFU/ml x 200ul	220	380	260
LVP004-SF	serum-free hSOX2 inducible particles	1x10E8 IFU/ml x 200ul	260	450	305
LVP005	hNANOG inducible lentiviral / lentivirus particles	5x10E7 IFU/ml x 200ul	220	380	260
LVP005-SF	serum-free hNANOG inducible particles	1x10E8 IFU/ml x 200ul	260	450	305
LVP006	hLIN28 inducible lentiviral / lentivirus particles	5x10E7 IFU/ml x 200ul	220	380	260
LVP006-SF	serum-free hLIN28 inducible particles	1x10E8 IFU/ml x 200ul	260	450	305
LVP007	hMyc inducible lentiviral / lentivirus particles	5x10E7 IFU/ml x 200ul	220	380	260
LVP007-SF	serum-free hMyc inducible particles	1x10E8 IFU/ml x 200ul	260	450	305
LVP008	hKLF4 inducible lentiviral / lentivirus particles	5x10E7 IFU/ml x 200ul	220	380	260
LVP008-SF	serum-free hKLF4 inducible particles	1x10E8 IFU/ml x 200ul	260	450	305

Human proteins (alphabetical)

Code	Description	Pack Size	GBP	CHF	EUR
LVP106	inducible expression lentiviral / lentivirus particles h aak1	1x10E7 IFU/ml x 200ul	290	510	350
LVP285	inducible expression lentiviral / lentivirus particles h ACTG1	1x10E7 IFU/ml x 200ul	290	510	350
LVP257	inducible expression lentiviral / lentivirus particles h acvr1	1x10E7 IFU/ml x 200ul	290	510	350
LVP118	inducible expression lentiviral / lentivirus particles h acvr1c	1x10E7 IFU/ml x 200ul	290	510	350
LVP116	inducible expression lentiviral / lentivirus particles h adck1	1x10E7 IFU/ml x 200ul	290	510	350
LVP031	inducible expression lentiviral / lentivirus particles h ak3	1x10E7 IFU/ml x 200ul	290	510	350
LVP117	inducible expression lentiviral / lentivirus particles h akt1	1x10E7 IFU/ml x 200ul	290	510	350
LVP121	inducible expression lentiviral / lentivirus particles h arrb2	1x10E7 IFU/ml x 200ul	290	510	350
LVP122	inducible expression lentiviral / lentivirus particles h atf2	1x10E7 IFU/ml x 200ul	290	510	350
LVP250	inducible expression lentiviral / lentivirus particles h aurka	1x10E7 IFU/ml x 200ul	290	510	350
LVP123	inducible expression lentiviral / lentivirus particles h aurkb	1x10E7 IFU/ml x 200ul	290	510	350
LVP124	inducible expression lentiviral / lentivirus particles h axl	1x10E7 IFU/ml x 200ul	290	510	350
LVP098	inducible expression lentiviral / lentivirus particles h b3galnt2	1x10E7 IFU/ml x 200ul	290	510	350
LVP288	inducible expression lentiviral / lentivirus particles h BAX	1x10E7 IFU/ml x 200ul	290	510	350
LVP126	inducible expression lentiviral / lentivirus particles h bckdk	1x10E7 IFU/ml x 200ul	290	510	350
LVP291	inducible expression lentiviral / lentivirus particles h BDKRB1	1x10E7 IFU/ml x 200ul	290	510	350
LVP036	inducible expression lentiviral / lentivirus particles h bmx	1x10E7 IFU/ml x 200ul	290	510	350
LVP078	inducible expression lentiviral / lentivirus particles h brd3	1x10E7 IFU/ml x 200ul	290	510	350
LVP247	inducible expression lentiviral / lentivirus particles h brsk2	1x10E7 IFU/ml x 200ul	290	510	350
LVP127	inducible expression lentiviral / lentivirus particles h btk	1x10E7 IFU/ml x 200ul	290	510	350
LVP043	inducible expression lentiviral / lentivirus particles h bub1	1x10E7 IFU/ml x 200ul	290	510	350
LVP129	inducible expression lentiviral / lentivirus particles h camk1	1x10E7 IFU/ml x 200ul	290	510	350
LVP130	inducible expression lentiviral / lentivirus particles h camk1d	1x10E7 IFU/ml x 200ul	290	510	350
LVP049	inducible expression lentiviral / lentivirus particles h camk2a	1x10E7 IFU/ml x 200ul	290	510	350

Code	Description	Pack Size	GBP	CHF	EUR
LVP133	inducible expression lentiviral / lentivirus particles h camk4	1x10E7 IFU/ml x 200ul	290	510	350
LVP203	inducible expression lentiviral / lentivirus particles h camkv	1x10E7 IFU/ml x 200ul	290	510	350
LVP261	inducible expression lentiviral / lentivirus particles h ccnd1	1x10E7 IFU/ml x 200ul	290	510	350
LVP134	inducible expression lentiviral / lentivirus particles h ccnd3	1x10E7 IFU/ml x 200ul	290	510	350
LVP135	inducible expression lentiviral / lentivirus particles h ccnh	1x10E7 IFU/ml x 200ul	290	510	350
LVP136	inducible expression lentiviral / lentivirus particles h cdc2l5	1x10E7 IFU/ml x 200ul	290	510	350
LVP137	inducible expression lentiviral / lentivirus particles h cdc42	1x10E7 IFU/ml x 200ul	290	510	350
LVP119	inducible expression lentiviral / lentivirus particles h cdk15	1x10E7 IFU/ml x 200ul	290	510	350
LVP232	inducible expression lentiviral / lentivirus particles h cdk16	1x10E7 IFU/ml x 200ul	290	510	350
LVP044	inducible expression lentiviral / lentivirus particles h cdk18	1x10E7 IFU/ml x 200ul	290	510	350
LVP263	inducible expression lentiviral / lentivirus particles h cdk2	1x10E7 IFU/ml x 200ul	290	510	350
LVP085	inducible expression lentiviral / lentivirus particles h cdk3	1x10E7 IFU/ml x 200ul	290	510	350
LVP264	inducible expression lentiviral / lentivirus particles h cdk4	1x10E7 IFU/ml x 200ul	290	510	350
LVP053	inducible expression lentiviral / lentivirus particles h cdk5r1	1x10E7 IFU/ml x 200ul	290	510	350
LVP265	inducible expression lentiviral / lentivirus particles h cdk7	1x10E7 IFU/ml x 200ul	290	510	350
LVP139	inducible expression lentiviral / lentivirus particles h cdk9	1x10E7 IFU/ml x 200ul	290	510	350
LVP079	inducible expression lentiviral / lentivirus particles h cdkl3	1x10E7 IFU/ml x 200ul	290	510	350
LVP082	inducible expression lentiviral / lentivirus particles h cdkl5	1x10E7 IFU/ml x 200ul	290	510	350
LVP140	inducible expression lentiviral / lentivirus particles h cdkn1a	1x10E7 IFU/ml x 200ul	290	510	350
LVP262	inducible expression lentiviral / lentivirus particles h cdkn2b	1x10E7 IFU/ml x 200ul	290	510	350
LVP142	inducible expression lentiviral / lentivirus particles h chek2	1x10E7 IFU/ml x 200ul	290	510	350
LVP113	inducible expression lentiviral / lentivirus particles h chuk	1x10E7 IFU/ml x 200ul	290	510	350
LVP035	inducible expression lentiviral / lentivirus particles h ckb	1x10E7 IFU/ml x 200ul	290	510	350
LVP055	inducible expression lentiviral / lentivirus particles h ckm	1x10E7 IFU/ml x 200ul	290	510	350
LVP065	inducible expression lentiviral / lentivirus particles h clcn2	1x10E7 IFU/ml x 200ul	290	510	350
LVP143	inducible expression lentiviral / lentivirus particles h clk1	1x10E7 IFU/ml x 200ul	290	510	350
LVP091	inducible expression lentiviral / lentivirus particles h clk3	1x10E7 IFU/ml x 200ul	290	510	350
LVP041	inducible expression lentiviral / lentivirus particles h cmpk1	1x10E7 IFU/ml x 200ul	290	510	350
LVP144	inducible expression lentiviral / lentivirus particles h col4a3bp	1x10E7 IFU/ml x 200ul	290	510	350
LVP145	inducible expression lentiviral / lentivirus particles h creb1	1x10E7 IFU/ml x 200ul	290	510	350
LVP200	inducible expression lentiviral / lentivirus particles h csnk1a1l	1x10E7 IFU/ml x 200ul	290	510	350
LVP146	inducible expression lentiviral / lentivirus particles h csnk1d	1x10E7 IFU/ml x 200ul	290	510	350
LVP038	inducible expression lentiviral / lentivirus particles h csnk1g1	1x10E7 IFU/ml x 200ul	290	510	350
LVP148	inducible expression lentiviral / lentivirus particles h csnk2a1	1x10E7 IFU/ml x 200ul	290	510	350
LVP282	inducible expression lentiviral / lentivirus particles h CTSD	1x10E7 IFU/ml x 200ul	290	510	350
LVP103	inducible expression lentiviral / lentivirus particles h cxcr4	1x10E7 IFU/ml x 200ul	290	510	350
LVP104	inducible expression lentiviral / lentivirus particles h cxcr7	1x10E7 IFU/ml x 200ul	290	510	350
LVP048	inducible expression lentiviral / lentivirus particles h dapk2	1x10E7 IFU/ml x 200ul	290	510	350
LVP061	inducible expression lentiviral / lentivirus particles h dclk1	1x10E7 IFU/ml x 200ul	290	510	350
LVP149	inducible expression lentiviral / lentivirus particles h ddit3	1x10E7 IFU/ml x 200ul	290	510	350
LVP150	inducible expression lentiviral / lentivirus particles h ddr1	1x10E7 IFU/ml x 200ul	290	510	350
LVP026	inducible expression lentiviral / lentivirus particles h dgke	1x10E7 IFU/ml x 200ul	290	510	350
LVP138	inducible expression lentiviral / lentivirus particles h dnajc28	1x10E7 IFU/ml x 200ul	290	510	350
LVP205	inducible expression lentiviral / lentivirus particles h dusp4	1x10E7 IFU/ml x 200ul	290	510	350
LVP206	inducible expression lentiviral / lentivirus particles h dusp6	1x10E7 IFU/ml x 200ul	290	510	350

Code	Description	Pack Size	GBP	CHF	EUR
LVP152	inducible expression lentiviral / lentivirus particles h dyrk2	1x10E7 IFU/ml x 200ul	290	510	350
LVP153	inducible expression lentiviral / lentivirus particles h eef2k	1x10E7 IFU/ml x 200ul	290	510	350
LVP071	inducible expression lentiviral / lentivirus particles h elavl1	1x10E7 IFU/ml x 200ul	290	510	350
LVP154	inducible expression lentiviral / lentivirus particles h elk1	1x10E7 IFU/ml x 200ul	290	510	350
LVP155	inducible expression lentiviral / lentivirus particles h elk4	1x10E7 IFU/ml x 200ul	290	510	350
LVP156	inducible expression lentiviral / lentivirus particles h epha2	1x10E7 IFU/ml x 200ul	290	510	350
LVP083	inducible expression lentiviral / lentivirus particles h epha4	1x10E7 IFU/ml x 200ul	290	510	350
LVP157	inducible expression lentiviral / lentivirus particles h ephb3	1x10E7 IFU/ml x 200ul	290	510	350
LVP267	inducible expression lentiviral / lentivirus particles h epo	1x10E7 IFU/ml x 200ul	290	510	350
LVP158	inducible expression lentiviral / lentivirus particles h erbb3	1x10E7 IFU/ml x 200ul	290	510	350
LVP159	inducible expression lentiviral / lentivirus particles h ets1	1x10E7 IFU/ml x 200ul	290	510	350
LVP160	inducible expression lentiviral / lentivirus particles h ets2	1x10E7 IFU/ml x 200ul	290	510	350
LVP161	inducible expression lentiviral / lentivirus particles h fastk	1x10E7 IFU/ml x 200ul	290	510	350
LVP296	inducible expression lentiviral / lentivirus particles h FGF18	1x10E7 IFU/ml x 200ul	290	510	350
LVP066	inducible expression lentiviral / lentivirus particles h fgfr1	1x10E7 IFU/ml x 200ul	290	510	350
LVP162	inducible expression lentiviral / lentivirus particles h fgfr2	1x10E7 IFU/ml x 200ul	290	510	350
LVP163	inducible expression lentiviral / lentivirus particles h fgr	1x10E7 IFU/ml x 200ul	290	510	350
LVP166	inducible expression lentiviral / lentivirus particles h fos	1x10E7 IFU/ml x 200ul	290	510	350
LVP258	inducible expression lentiviral / lentivirus particles h foxo3	1x10E7 IFU/ml x 200ul	290	510	350
LVP077	inducible expression lentiviral / lentivirus particles h fyn	1x10E7 IFU/ml x 200ul	290	510	350
LVP039	inducible expression lentiviral / lentivirus particles h gak	1x10E7 IFU/ml x 200ul	290	510	350
LVP040	inducible expression lentiviral / lentivirus particles h galk2	1x10E7 IFU/ml x 200ul	290	510	350
LVP074	inducible expression lentiviral / lentivirus particles h gba	1x10E7 IFU/ml x 200ul	290	510	350
LVP045	inducible expression lentiviral / lentivirus particles h gck	1x10E7 IFU/ml x 200ul	290	510	350
LVP266	inducible expression lentiviral / lentivirus particles h gdf9	1x10E7 IFU/ml x 200ul	290	510	350
LVP273	inducible expression lentiviral / lentivirus particles h GGA3	1x10E7 IFU/ml x 200ul	290	510	350
LVP015	inducible expression lentiviral / lentivirus particles h glul	1x10E7 IFU/ml x 200ul	290	510	350
LVP076	inducible expression lentiviral / lentivirus particles h gorasp2	1x10E7 IFU/ml x 200ul	290	510	350
LVP112	inducible expression lentiviral / lentivirus particles h gpsm3	1x10E7 IFU/ml x 200ul	290	510	350
LVP046	inducible expression lentiviral / lentivirus particles h grk6	1x10E7 IFU/ml x 200ul	290	510	350
LVP167	inducible expression lentiviral / lentivirus particles h gsk3b	1x10E7 IFU/ml x 200ul	290	510	350
LVP277	inducible expression lentiviral / lentivirus particles h IFNG	1x10E7 IFU/ml x 200ul	290	510	350
LVP169	inducible expression lentiviral / lentivirus particles h ikbkb	1x10E7 IFU/ml x 200ul	290	510	350
LVP170	inducible expression lentiviral / lentivirus particles h ikbkg	1x10E7 IFU/ml x 200ul	290	510	350
LVP073	inducible expression lentiviral / lentivirus particles h il17br	1x10E7 IFU/ml x 200ul	290	510	350
LVP293	inducible expression lentiviral / lentivirus particles h IL2	1x10E7 IFU/ml x 200ul	290	510	350
LVP295	inducible expression lentiviral / lentivirus particles h IL20	1x10E7 IFU/ml x 200ul	290	510	350
LVP292	inducible expression lentiviral / lentivirus particles h IL3	1x10E7 IFU/ml x 200ul	290	510	350
LVP209	inducible expression lentiviral / lentivirus particles h irak3	1x10E7 IFU/ml x 200ul	290	510	350
LVP172	inducible expression lentiviral / lentivirus particles h jun	1x10E7 IFU/ml x 200ul	290	510	350
LVP290	inducible expression lentiviral / lentivirus particles h KCNJ5	1x10E7 IFU/ml x 200ul	290	510	350
LVP294	inducible expression lentiviral / lentivirus particles h KITLG	1x10E7 IFU/ml x 200ul	290	510	350
LVP101	inducible expression lentiviral / lentivirus particles h lbr	1x10E7 IFU/ml x 200ul	290	510	350
LVP033	inducible expression lentiviral / lentivirus particles h limk2	1x10E7 IFU/ml x 200ul	290	510	350
LVP174	inducible expression lentiviral / lentivirus particles h lmna	1x10E7 IFU/ml x 200ul	290	510	350

Code	Description	Pack Size	GBP	CHF	EUR
LVP092	inducible expression lentiviral / lentivirus particles h mak	1x10E7 IFU/ml x 200ul	290	510	350
LVP176	inducible expression lentiviral / lentivirus particles h map2k2	1x10E7 IFU/ml x 200ul	290	510	350
LVP177	inducible expression lentiviral / lentivirus particles h map2k3	1x10E7 IFU/ml x 200ul	290	510	350
LVP178	inducible expression lentiviral / lentivirus particles h map2k4	1x10E7 IFU/ml x 200ul	290	510	350
LVP228	inducible expression lentiviral / lentivirus particles h map2k5	1x10E7 IFU/ml x 200ul	290	510	350
LVP180	inducible expression lentiviral / lentivirus particles h map3k12	1x10E7 IFU/ml x 200ul	290	510	350
LVP181	inducible expression lentiviral / lentivirus particles h map3k14	1x10E7 IFU/ml x 200ul	290	510	350
LVP182	inducible expression lentiviral / lentivirus particles h map3k5	1x10E7 IFU/ml x 200ul	290	510	350
LVP183	inducible expression lentiviral / lentivirus particles h map3k7	1x10E7 IFU/ml x 200ul	290	510	350
LVP184	inducible expression lentiviral / lentivirus particles h map3k8	1x10E7 IFU/ml x 200ul	290	510	350
LVP185	inducible expression lentiviral / lentivirus particles h mapk1	1x10E7 IFU/ml x 200ul	290	510	350
LVP186	inducible expression lentiviral / lentivirus particles h mapk12	1x10E7 IFU/ml x 200ul	290	510	350
LVP187	inducible expression lentiviral / lentivirus particles h mapk13	1x10E7 IFU/ml x 200ul	290	510	350
LVP188	inducible expression lentiviral / lentivirus particles h mapk3	1x10E7 IFU/ml x 200ul	290	510	350
LVP189	inducible expression lentiviral / lentivirus particles h mapk6	1x10E7 IFU/ml x 200ul	290	510	350
LVP190	inducible expression lentiviral / lentivirus particles h mapk7	1x10E7 IFU/ml x 200ul	290	510	350
LVP191	inducible expression lentiviral / lentivirus particles h mapk8	1x10E7 IFU/ml x 200ul	290	510	350
LVP192	inducible expression lentiviral / lentivirus particles h mapk9	1x10E7 IFU/ml x 200ul	290	510	350
LVP075	inducible expression lentiviral / lentivirus particles h mapkapk3	1x10E7 IFU/ml x 200ul	290	510	350
LVP195	inducible expression lentiviral / lentivirus particles h mapkapk5	1x10E7 IFU/ml x 200ul	290	510	350
LVP032	inducible expression lentiviral / lentivirus particles h mark2	1x10E7 IFU/ml x 200ul	290	510	350
LVP029	inducible expression lentiviral / lentivirus particles h mark3	1x10E7 IFU/ml x 200ul	290	510	350
LVP196	inducible expression lentiviral / lentivirus particles h mast2	1x10E7 IFU/ml x 200ul	290	510	350
LVP030	inducible expression lentiviral / lentivirus particles h matk	1x10E7 IFU/ml x 200ul	290	510	350
LVP197	inducible expression lentiviral / lentivirus particles h max	1x10E7 IFU/ml x 200ul	290	510	350
LVP198	inducible expression lentiviral / lentivirus particles h mef2b	1x10E7 IFU/ml x 200ul	290	510	350
LVP199	inducible expression lentiviral / lentivirus particles h mef2c	1x10E7 IFU/ml x 200ul	290	510	350
LVP204	inducible expression lentiviral / lentivirus particles h mknk1	1x10E7 IFU/ml x 200ul	290	510	350
LVP229	inducible expression lentiviral / lentivirus particles h mre11a	1x10E7 IFU/ml x 200ul	290	510	350
LVP207	inducible expression lentiviral / lentivirus particles h myd88	1x10E7 IFU/ml x 200ul	290	510	350
LVP111	inducible expression lentiviral / lentivirus particles h nbn	1x10E7 IFU/ml x 200ul	290	510	350
LVP286	inducible expression lentiviral / lentivirus particles h NCOA4	1x10E7 IFU/ml x 200ul	290	510	350
LVP165	inducible expression lentiviral / lentivirus particles h nek10	1x10E7 IFU/ml x 200ul	290	510	350
LVP089	inducible expression lentiviral / lentivirus particles h nek11	1x10E7 IFU/ml x 200ul	290	510	350
LVP230	inducible expression lentiviral / lentivirus particles h nek6	1x10E7 IFU/ml x 200ul	290	510	350
LVP208	inducible expression lentiviral / lentivirus particles h nfatc3	1x10E7 IFU/ml x 200ul	290	510	350
LVP100	inducible expression lentiviral / lentivirus particles h nfic	1x10E7 IFU/ml x 200ul	290	510	350
LVP210	inducible expression lentiviral / lentivirus particles h nfkb1	1x10E7 IFU/ml x 200ul	290	510	350
LVP168	inducible expression lentiviral / lentivirus particles h nfkbia	1x10E7 IFU/ml x 200ul	290	510	350
LVP201	inducible expression lentiviral / lentivirus particles h nim1	1x10E7 IFU/ml x 200ul	290	510	350
LVP050	inducible expression lentiviral / lentivirus particles h nme1	1x10E7 IFU/ml x 200ul	290	510	350
LVP028	inducible expression lentiviral / lentivirus particles h nme5	1x10E7 IFU/ml x 200ul	290	510	350
LVP110	inducible expression lentiviral / lentivirus particles h npm1	1x10E7 IFU/ml x 200ul	290	510	350
LVP108	inducible expression lentiviral / lentivirus particles h nr3c1 variant 1	1x10E7 IFU/ml x 200ul	290	510	350
LVP109	inducible expression lentiviral / lentivirus particles h nr3c1 variant 6	1x10E7 IFU/ml x 200ul	290	510	350

Code	Description	Pack Size	GBP	CHF	EUR
LVP272	inducible expression lentiviral / lentivirus particles h NRAS	1x10E7 IFU/ml x 200ul	290	510	350
LVP047	inducible expression lentiviral / lentivirus particles h nuak2	1x10E7 IFU/ml x 200ul	290	510	350
LVP259	inducible expression lentiviral / lentivirus particles h oxsr1	1x10E7 IFU/ml x 200ul	290	510	350
LVP211	inducible expression lentiviral / lentivirus particles h pak2	1x10E7 IFU/ml x 200ul	290	510	350
LVP231	inducible expression lentiviral / lentivirus particles h pak4	1x10E7 IFU/ml x 200ul	290	510	350
LVP037	inducible expression lentiviral / lentivirus particles h pak6	1x10E7 IFU/ml x 200ul	290	510	350
LVP125	inducible expression lentiviral / lentivirus particles h pars2	1x10E7 IFU/ml x 200ul	290	510	350
LVP095	inducible expression lentiviral / lentivirus particles h pbk	1x10E7 IFU/ml x 200ul	290	510	350
LVP278	inducible expression lentiviral / lentivirus particles h PCBP2	1x10E7 IFU/ml x 200ul	290	510	350
LVP102	inducible expression lentiviral / lentivirus particles h pdha1	1x10E7 IFU/ml x 200ul	290	510	350
LVP052	inducible expression lentiviral / lentivirus particles h pdik1l	1x10E7 IFU/ml x 200ul	290	510	350
LVP060	inducible expression lentiviral / lentivirus particles h pdk3	1x10E7 IFU/ml x 200ul	290	510	350
LVP212	inducible expression lentiviral / lentivirus particles h pdk4	1x10E7 IFU/ml x 200ul	290	510	350
LVP025	inducible expression lentiviral / lentivirus particles h pfkml	1x10E7 IFU/ml x 200ul	290	510	350
LVP213	inducible expression lentiviral / lentivirus particles h phkg2	1x10E7 IFU/ml x 200ul	290	510	350
LVP214	inducible expression lentiviral / lentivirus particles h pik3cb	1x10E7 IFU/ml x 200ul	290	510	350
LVP215	inducible expression lentiviral / lentivirus particles h pik3r3	1x10E7 IFU/ml x 200ul	290	510	350
LVP062	inducible expression lentiviral / lentivirus particles h pim1	1x10E7 IFU/ml x 200ul	290	510	350
LVP034	inducible expression lentiviral / lentivirus particles h pim2	1x10E7 IFU/ml x 200ul	290	510	350
LVP058	inducible expression lentiviral / lentivirus particles h pip5k1a	1x10E7 IFU/ml x 200ul	290	510	350
LVP233	inducible expression lentiviral / lentivirus particles h pkmyt1	1x10E7 IFU/ml x 200ul	290	510	350
LVP234	inducible expression lentiviral / lentivirus particles h plk1	1x10E7 IFU/ml x 200ul	290	510	350
LVP226	inducible expression lentiviral / lentivirus particles h plk4	1x10E7 IFU/ml x 200ul	290	510	350
LVP068	inducible expression lentiviral / lentivirus particles h pomp	1x10E7 IFU/ml x 200ul	290	510	350
LVP216	inducible expression lentiviral / lentivirus particles h pparg	1x10E7 IFU/ml x 200ul	290	510	350
LVP072	inducible expression lentiviral / lentivirus particles h ppp1r7	1x10E7 IFU/ml x 200ul	290	510	350
LVP217	inducible expression lentiviral / lentivirus particles h prkaa1	1x10E7 IFU/ml x 200ul	290	510	350
LVP094	inducible expression lentiviral / lentivirus particles h prkaa2	1x10E7 IFU/ml x 200ul	290	510	350
LVP218	inducible expression lentiviral / lentivirus particles h prkaca	1x10E7 IFU/ml x 200ul	290	510	350
LVP235	inducible expression lentiviral / lentivirus particles h prkacb	1x10E7 IFU/ml x 200ul	290	510	350
LVP219	inducible expression lentiviral / lentivirus particles h prkca	1x10E7 IFU/ml x 200ul	290	510	350
LVP236	inducible expression lentiviral / lentivirus particles h prkcg	1x10E7 IFU/ml x 200ul	290	510	350
LVP099	inducible expression lentiviral / lentivirus particles h prom1	1x10E7 IFU/ml x 200ul	290	510	350
LVP220	inducible expression lentiviral / lentivirus particles h ptk2b	1x10E7 IFU/ml x 200ul	290	510	350
LVP088	inducible expression lentiviral / lentivirus particles h rage	1x10E7 IFU/ml x 200ul	290	510	350
LVP274	inducible expression lentiviral / lentivirus particles h RaBP1	1x10E7 IFU/ml x 200ul	290	510	350
LVP275	inducible expression lentiviral / lentivirus particles h RanBP1	1x10E7 IFU/ml x 200ul	290	510	350
LVP284	inducible expression lentiviral / lentivirus particles h RASSF5	1x10E7 IFU/ml x 200ul	290	510	350
LVP237	inducible expression lentiviral / lentivirus particles h rela	1x10E7 IFU/ml x 200ul	290	510	350
LVP221	inducible expression lentiviral / lentivirus particles h ret	1x10E7 IFU/ml x 200ul	290	510	350
LVP280	inducible expression lentiviral / lentivirus particles h RGS16	1x10E7 IFU/ml x 200ul	290	510	350
LVP251	inducible expression lentiviral / lentivirus particles h riok3	1x10E7 IFU/ml x 200ul	290	510	350
LVP222	inducible expression lentiviral / lentivirus particles h ripk2	1x10E7 IFU/ml x 200ul	290	510	350
LVP114	inducible expression lentiviral / lentivirus particles h rnaseh2a	1x10E7 IFU/ml x 200ul	290	510	350
LVP064	inducible expression lentiviral / lentivirus particles h rp2	1x10E7 IFU/ml x 200ul	290	510	350

Code	Description	Pack Size	GBP	CHF	EUR
LVP193	inducible expression lentiviral / lentivirus particles h rps6ka1	1x10E7 IFU/ml x 200ul	290	510	350
LVP238	inducible expression lentiviral / lentivirus particles h rps6ka2	1x10E7 IFU/ml x 200ul	290	510	350
LVP194	inducible expression lentiviral / lentivirus particles h rps6ka3	1x10E7 IFU/ml x 200ul	290	510	350
LVP239	inducible expression lentiviral / lentivirus particles h rps6kb1	1x10E7 IFU/ml x 200ul	290	510	350
LVP086	inducible expression lentiviral / lentivirus particles h rps6kl1	1x10E7 IFU/ml x 200ul	290	510	350
LVP281	inducible expression lentiviral / lentivirus particles h RUVBL2	1x10E7 IFU/ml x 200ul	290	510	350
LVP223	inducible expression lentiviral / lentivirus particles h sgk2	1x10E7 IFU/ml x 200ul	290	510	350
LVP240	inducible expression lentiviral / lentivirus particles h sgk3	1x10E7 IFU/ml x 200ul	290	510	350
LVP175	inducible expression lentiviral / lentivirus particles h smad2	1x10E7 IFU/ml x 200ul	290	510	350
LVP241	inducible expression lentiviral / lentivirus particles h sp1	1x10E7 IFU/ml x 200ul	290	510	350
LVP242	inducible expression lentiviral / lentivirus particles h spag9	1x10E7 IFU/ml x 200ul	290	510	350
LVP120	inducible expression lentiviral / lentivirus particles h speg	1x10E7 IFU/ml x 200ul	290	510	350
LVP132	inducible expression lentiviral / lentivirus particles h spz1	1x10E7 IFU/ml x 200ul	290	510	350
LVP243	inducible expression lentiviral / lentivirus particles h srp1	1x10E7 IFU/ml x 200ul	290	510	350
LVP224	inducible expression lentiviral / lentivirus particles h stat3	1x10E7 IFU/ml x 200ul	290	510	350
LVP081	inducible expression lentiviral / lentivirus particles h stk10	1x10E7 IFU/ml x 200ul	290	510	350
LVP090	inducible expression lentiviral / lentivirus particles h stk11	1x10E7 IFU/ml x 200ul	290	510	350
LVP225	inducible expression lentiviral / lentivirus particles h stk16	1x10E7 IFU/ml x 200ul	290	510	350
LVP056	inducible expression lentiviral / lentivirus particles h stk17b	1x10E7 IFU/ml x 200ul	290	510	350
LVP246	inducible expression lentiviral / lentivirus particles h stk25	1x10E7 IFU/ml x 200ul	290	510	350
LVP248	inducible expression lentiviral / lentivirus particles h stk3	1x10E7 IFU/ml x 200ul	290	510	350
LVP063	inducible expression lentiviral / lentivirus particles h stk32c	1x10E7 IFU/ml x 200ul	290	510	350
LVP249	inducible expression lentiviral / lentivirus particles h stk38	1x10E7 IFU/ml x 200ul	290	510	350
LVP080	inducible expression lentiviral / lentivirus particles h stk4	1x10E7 IFU/ml x 200ul	290	510	350
LVP202	inducible expression lentiviral / lentivirus particles h stk40	1x10E7 IFU/ml x 200ul	290	510	350
LVP107	inducible expression lentiviral / lentivirus particles h stradb	1x10E7 IFU/ml x 200ul	290	510	350
LVP260	inducible expression lentiviral / lentivirus particles h stradb	1x10E7 IFU/ml x 200ul	290	510	350
LVP151	inducible expression lentiviral / lentivirus particles h styk1	1x10E7 IFU/ml x 200ul	290	510	350
LVP171	inducible expression lentiviral / lentivirus particles h taok3	1x10E7 IFU/ml x 200ul	290	510	350
LVP276	inducible expression lentiviral / lentivirus particles h TAPBP	1x10E7 IFU/ml x 200ul	290	510	350
LVP105	inducible expression lentiviral / lentivirus particles h tert	1x10E7 IFU/ml x 200ul	290	510	350
LVP289	inducible expression lentiviral / lentivirus particles h TGFB1	1x10E7 IFU/ml x 200ul	290	510	350
LVP054	inducible expression lentiviral / lentivirus particles h tgfr2	1x10E7 IFU/ml x 200ul	290	510	350
LVP097	inducible expression lentiviral / lentivirus particles h tmem50b	1x10E7 IFU/ml x 200ul	290	510	350
LVP287	inducible expression lentiviral / lentivirus particles h TNF	1x10E7 IFU/ml x 200ul	290	510	350
LVP173	inducible expression lentiviral / lentivirus particles h tnk	1x10E7 IFU/ml x 200ul	290	510	350
LVP115	inducible expression lentiviral / lentivirus particles h tnk2	1x10E7 IFU/ml x 200ul	290	510	350
LVP253	inducible expression lentiviral / lentivirus particles h tp53	1x10E7 IFU/ml x 200ul	290	510	350
LVP069	inducible expression lentiviral / lentivirus particles h tpt1	1x10E7 IFU/ml x 200ul	290	510	350
LVP244	inducible expression lentiviral / lentivirus particles h tpte	1x10E7 IFU/ml x 200ul	290	510	350
LVP042	inducible expression lentiviral / lentivirus particles h trib2	1x10E7 IFU/ml x 200ul	290	510	350
LVP128	inducible expression lentiviral / lentivirus particles h trib3	1x10E7 IFU/ml x 200ul	290	510	350
LVP084	inducible expression lentiviral / lentivirus particles h tssk1b	1x10E7 IFU/ml x 200ul	290	510	350
LVP227	inducible expression lentiviral / lentivirus particles h tssk2	1x10E7 IFU/ml x 200ul	290	510	350
LVP245	inducible expression lentiviral / lentivirus particles h tssk3	1x10E7 IFU/ml x 200ul	290	510	350

Code	Description	Pack Size	GBP	CHF	EUR
LVP254	inducible expression lentiviral / lentivirus particles h ttbk2	1x10E7 IFU/ml x 200ul	290	510	350
LVP096	inducible expression lentiviral / lentivirus particles h tuba1b	1x10E7 IFU/ml x 200ul	290	510	350
LVP087	inducible expression lentiviral / lentivirus particles h twf1	1x10E7 IFU/ml x 200ul	290	510	350
LVP093	inducible expression lentiviral / lentivirus particles h twf2	1x10E7 IFU/ml x 200ul	290	510	350
LVP067	inducible expression lentiviral / lentivirus particles h txn	1x10E7 IFU/ml x 200ul	290	510	350
LVP051	inducible expression lentiviral / lentivirus particles h tyk2	1x10E7 IFU/ml x 200ul	290	510	350
LVP022	inducible expression lentiviral / lentivirus particles h ubd	1x10E7 IFU/ml x 200ul	290	510	350
LVP164	inducible expression lentiviral / lentivirus particles h ulk4	1x10E7 IFU/ml x 200ul	290	510	350
LVP255	inducible expression lentiviral / lentivirus particles h vrk3	1x10E7 IFU/ml x 200ul	290	510	350
LVP252	inducible expression lentiviral / lentivirus particles h wee1	1x10E7 IFU/ml x 200ul	290	510	350
LVP057	inducible expression lentiviral / lentivirus particles h wif1	1x10E7 IFU/ml x 200ul	290	510	350
LVP070	inducible expression lentiviral / lentivirus particles h xpo5	1x10E7 IFU/ml x 200ul	290	510	350
LVP059	inducible expression lentiviral / lentivirus particles h yes1	1x10E7 IFU/ml x 200ul	290	510	350
LVP256	inducible expression lentiviral / lentivirus particles h zap70	1x10E7 IFU/ml x 200ul	290	510	350
LVP279	inducible expression lentiviral / lentivirus particles h ZFP36	1x10E7 IFU/ml x 200ul	290	510	350

Mouse proteins for inducing pluripotent stem cells (iPSC)

Code	Description	Pack Size	GBP	CHF	EUR
LVP-stems-m	Full set (6 mouse stem factors) serum-free particles	200ul x 6	1205	2100	1435
LVP003m	m oct4 inducible lentiviral / lentivirus particles (mouse)	1x10E7 IFU/ml x 200ul	260	450	305
LVP004m	m sox2 inducible lentiviral / lentivirus particles (mouse)	1x10E7 IFU/ml x 200ul	260	450	305
LVP005m	m nanog inducible lentiviral / lentivirus particles (mouse)	5x10E7 IFU/ml x 200ul	220	380	260
LVP006m	m lin28 inducible lentiviral / lentivirus particles (mouse)	5x10E7 IFU/ml x 200ul	220	380	260
LVP007m	m myc inducible lentiviral / lentivirus particles (mouse)	5x10E7 IFU/ml x 200ul	220	380	260
LVP008m	m klf4 inducible lentiviral / lentivirus particles (mouse)	5x10E7 IFU/ml x 200ul	220	380	260

Mouse proteins (alphabetical)

Code	Description	Pack Size	GBP	CHF	EUR
LVP131	m acp2 inducible lentiviral / lentivirus particles (mouse)	1x10E7 IFU/ml x 200ul	290	510	350
LVP179	m csn3 inducible lentiviral / lentivirus particles (mouse)	1x10E7 IFU/ml x 200ul	290	510	350
LVP141	m gadd45a inducible lentiviral / lentivirus particles (mouse)	1x10E7 IFU/ml x 200ul	290	510	350
LVP270	m mef2c inducible lentiviral / lentivirus particles (mouse)	1x10E7 IFU/ml x 200ul	290	510	350
LVP271	m mesp1 inducible lentiviral / lentivirus particles (mouse)	1x10E7 IFU/ml x 200ul	290	510	350
LVP268	m mitf inducible lentiviral / lentivirus particles (mouse)	1x10E7 IFU/ml x 200ul	290	510	350
LVP147	m myd88 inducible lentiviral / lentivirus particles (mouse)	1x10E7 IFU/ml x 200ul	290	510	350
LVP269	m pax6 inducible lentiviral / lentivirus particles (mouse)	1x10E7 IFU/ml x 200ul	290	510	350



UK & Rest of World

184 Milton Park,
Abingdon, OX14 4SE - UK
Tel: +44 (0) 1235 828 200
Fax: +44 (0) 1235 820 482

Switzerland

Centro Nord-Sud 2E
CH-6934 Bioggio-Lugano
Tel: +41 (0) 91 604 55 22
Fax: +41 (0) 91 605 17 85

Deutschland

Tel: +49 (0) 69 779099
Fax: +49 (0) 69 13376880

amsbio

info@amsbio.com

AMS Biotechnology (Europe) Ltd, Registered
in England & Wales No.2117791