

PURIFIED PROTEINS

PROTEIN : HUMAN IgG1 kappa

CAT. NO : 523

LOT NO : 2400 EXP: DEC. '13

SOURCE : IgG1 kappa Myeloma serum *

PURITY : 98+%

SIZE : 1.0 mg

BUFFER : 1.0 ml 0.15M NaCl
0.01M phosphate
0.3 M glycine
0.1% sodium azide pH 7.2

PROT. CONC. : 1.0 mg/ml

STORAGE : +4°C for continuous use
-20°C for extended storage
Avoid repeated freezing and thawing,
aliquot if necessary

* Shown to be negative for HBsAg and antibodies to HCV and HIV-1/2 by a certified test.

This product is for laboratory use only.
Not for human use.



UK & Rest of World
184 Milton Park, Abingdon
OX14 4SE, Oxon, UK
Tel: +44 (0) 1235 828 200
Fax: +44 (0) 1235 820 482

Switzerland
Centro Nord-Sud 2E
CH-6934 Bioggio-Lugano
Tel: +41 (0) 91 604 55 22
Fax: +41 (0) 91 605 17 85

Deutschland
Bockenheimer Landstr. 17/19
60325 Frankfurt/Main
Tel: +49 (0) 69 779099
Fax: +49 (0) 69 13376880

North America
23591 El Toro Rd, Suite #167
Lake Forest, CA 92630
Tel: + 1 800 987 0985
Fax: + 1 949 265 7703