

## PURIFIED PROTEINS

PROTEIN : HUMAN IgM kappa

CAT. NO : 541

LOT NO : 2426 EXP: DEC. '13

SOURCE : IgM kappa myeloma serum \*

PURITY : 98+%

SIZE : 1.0 mg


BUFFER : 1.0 ml 0.15M NaCl  
0.01M phosphate  
1.0 M glycine  
0.1% sodium azide pH 7.2

PROT. CONC. : 1.0 mg/ml

STORAGE : +4°C, DO NOT FREEZE

\* Shown to be negative for HBsAg and antibodies to HCV and HIV-1/2 by a certified test.

This product is for laboratory use only.  
Not for human use.



**United States**  
23591 El Toro Rd, Suite #167  
Lake Forest, CA 92630  
Toll Free: + 1 800 987 0985  
Tel: + 1 949 265 7717  
Fax: + 1 949 265 7703

**UK & Rest of World**  
184 Milton Park, Abingdon  
OX14 4SE, Oxon, UK  
Tel: +44 (0) 1235 828 200  
Fax: +44 (0) 1235 820 482

**Switzerland**  
Centro Nord-Sud 2E  
CH-6934 Bioggio-Lugano  
Tel: +41 (0) 91 604 55 22  
Fax: +41 (0) 91 605 17 85

**Deutschland**  
Bockenheimer Landstr. 17/19  
60325 Frankfurt/Main  
Tel: +49 (0) 69 779099  
Fax: +49 (0) 69 13376880