

Human primary tumour and adjacent normal paraffin tissue sections

Storage Conditions: Cool dry at room temperature. Keep at 4°C for long term storage. One year from the date of receipt under proper storage condition.

Applications: Human matched pair paraffin tissue sections are ideal for:

- Cellular localization of specific mRNA and protein expression in tumour tissue
- Identification of tumour specific genes, proteins by *in situ* hybridization, *in situ* PCR detection and immunohistochemistry
- Comparison of the locations of novel tumour genes in different human tissues.

Quality Control: At least one of the tissue slides from each lot was stained with H & E to ensure the quality.

Description: The human tumour and adjacent sections (PP) consist of paraffin samples from the same donor. The tissues were fixed in formalin immediately after excision and embedded in paraffin. The tissue section with 5-10 µm thickness is mounted on positively charged glass slides and fixed by cold acetone. 5 slides of primary tumour and another 5 slides of adjacent normal tissue are provided.

Shipping: Room temperature.

Intended use of the products: All products sold by AMSBIO are intended for research use only and must not be used for medical diagnostic or drug purposes. AMSBIO products may not resell, modified for resale or used to manufacture commercial products without written approval from AMSBIO.

Limited product warranty: AMSBIO has a total quality assurance policy. If any product does not meet the specification stated in the product data sheet will be replaced at no charge. This warranty limits AMSBIO liability to the replacement of the product only. No other warranties of any kind, express or implied, including without limitation, implied warranties of merchantability or fitness for a particular purpose, are provided by AMSBIO. AMSBIO shall not be liable for any incidental, direct, special or consequential damages, even if it is aware of the possibility of such damages.