

ImmunoHistoChemistry kit

Quantity: 500 slides

Applications: This kit can be used to detect cell surface and intracellular antigens in frozen or paraffin-embedded tissue sections mounted on slides, floating sections in culture plates and immobilized cells with high sensitivity and low background.

Description: Kits are easy to use, provided with instructions for each step from preparation of tissue sections to mounting slides with coverslips. In addition to the immunostaining procedure, the user manual also contains detailed technical information about processing frozen and paraffin tissue sections for immune-staining, antigen retrieval and trouble-shootings.

Content: The kit contains reagents required for immunostaining of 500 slides:

- Serum blocking buffer
- Antibody diluent
- Biotinylated secondary antibody
- Hematoxylin
- Mounting medium
- Streptavidin-peroxidase conjugate*
- Streptavidin-alkaline phosphatase conjugate*
- Streptavidin-FITC conjugate*
- DAB substrate*
- Peroxidase suppresser*
- Alkaline phosphatase suppresser*
- Substrate*

* This reagents are provided depending on the immunohistochemistry method.

Intended use of the products: All products sold by AMSBIO are intended for research use only and must not be used for medical diagnostic or drug purposes. AMSBIO products may not resell, modified for resale or used to manufacture commercial products without written approval from AMSBIO.

Limited product warranty: AMSBIO has a total quality assurance policy. If an y product does not meet the specification stated in the product data sheet will be replaced at no charge. This warranty limits AMSBIO liability to the replacement of the product only. No other warranties of an y kind, express or implied, including without limitation, implied warranties of merchantability or fitness for a particular purpose, are provided by AMSBIO. AMSBIO shall not liable for an y incidental, in direct, special or consequential damages, even if it is aware of the possibility of such damages.